

Advanced Elderly Fully-Functional Nursing Manikin (Male) Supplied by Desco Medical India



Specification:-

Multi-functional Old Female lifelike Manikin.
Specifically Designed For Professional Nursing Training.
High Quality Manikin Provides Wide Range Of Exercise
Possible Including Both First And Emergency Situation.

NURSING SKILL:

1. Hair and Face Washing
2. Pupil Observation, Bilateral Pupil States: Normal,
3. Eye and Ear Washing and Administering
4. Removal and Insertion of Hearing Aid
5. Oral Cavity and Artificial Teeth Care
6. Endotracheal Intubation, Automatically Monitor
7. Incision Of Trachea Nursing
8. Sputum Suction
9. Oxygen Inhaling
10. Oral and Nasal Feeding

11. Gastrolavage
12. Venipuncture, Injection, Transfusion (Arm)
13. Deltoid Subcutaneous Injection Steps Of 1 MM Hg.
14. Thoracic Cavity, Abdominal Cavity, Liver, Bone, Marrow and Lumbar Puncture
15. Breast Nursing, Mammary Gland Examination
16. Enema people while nursing.
17. Female Urethral Catheterization
18. Female Bladder Irrigation
19. Stoma Drainage
20. Colostomy
- 21 . Buttocks Intramuscular Injection
22. Large Area Of Sacral Ulcer
23. Cancer Tumor Contrast
24. Cleavage Lines Contrast Mydrisis
25. Holistic Nursing: Sponge Bath, Replacing Clothes Cold and Hot Therapy
26. Realistic Limbs Joint, Stiff Simulative Joints
27. Truck May Lean Forward, Can Sit It The Wheelchair
 - a. Truck-Rotation, Flexibility Weather The Tracheal
 - b. Neck-Rotation, Flexibility, Lateral e rvature Intubation Is Correct Or Not Through Sensor, Support
 - c. Shoulder and Buttock-Endoduction, Abduction, Rotation, Auscultation to Detect Intubation Site. Flexibility
 - d. Elbow-Internal Rotation, External Rotation,
 - e. Knee-Internal Rotation, External Rotation
 - f. Wrist-Internal Rotation, Flexibility, Stretching, Bending
 - g. Ankle-Introversion, Eversion Endoduction, Abduction
28. Non-Invasive Blood Pressure Measurement
 - a. Systolic and Diastolic Pressure May Be Set Individually In
 - b. Systolic Pressure 0-300MM Hg, Diastolic Pressure 0-300Hg
 - c. Korotkoff Sound Volume In Ten Steps, 0-9

This Manikin helps nursing students to understand old age

We are Exporters and Suppliers of Advanced Elderly Fully-Functional Nursing Manikin (Male) at wholesale competitive prices from India.

MODEL NO:
MMNT 485

Advanced Elderly Fully-Functional Nursing Manikin (Male)