



Autoclave Pressure Cooker Powder Coated Non Electric (MODEL NO : AUPC 251)



Autoclave Pressure Cooker Powder Coated Non Electric Manufactured by Desco Medical India



Decide with Confidence

Autoclave Pressure Cooker Powder Coated Non
Electric



www.descomedicalindia.com

www.autoclavesterilizersindia.com

www.descoonline.in

Autoclave Pressure Cooker Powder Coated Non Electric Manufactured by Desco Medical India

Specifications :-

Autoclave Pressure Cooker is also known as Pressure Steam Sterilizers, Steam Autoclave, Portable Autoclaves Sterilizers
 Autoclave Pressure Cooker is a part of Steam Sterilizing Equipment
 Autoclave Pressure Cooker is Non Electric and is Epoxy Powder Coated.
 The body and lid has Seamless Construction from Extra Thick Aluminium Sheet Body & Lid for more safety and life long performance.
 The Anti-Bulging & Fuel Saver Base is provided for Additional Strength, Fuel Saving and Time Saving.
 The depression at centre of the lid helps in Additional Strength.
 The Gasket is ISI Marked food grade synthetic rubber gasket for better reliability, hygiene & longer life.
 The Gasket Release System (GRS) & fusible plug type safety valve is available for Triple-Safety.
 The pressure gauge is Colour-coded Pressure Gauge to indicate the safe and danger zones.
 The repeating type safety valve is provided for better Safety and longer & error-free performance.
 The vacuum release & steam release valves are made of Stainless Steel for longer life and to avoid toxic deposits on valves due to oxidation.
 The handles is made from heat resistant, impact tolerant, finest quality bakelite, fitted with Stainless Steel screws to avoid easy loosening of handles.
 The inside tripod stand is made of Stainless Steel with extra height to avoid wetting of instruments.
 The outer stand is well designed to provide proper stability to the Autoclave Pressure Cooker.
 We are Manufacturer, Exporters and Suppliers of Autoclave Pressure Cooker Powder Coated Non Electric at wholesale competitive prices from India.

Model No	Size in inch	Size in mm	Capacity (Ltr)
AUPC 251	9 x 9 Inch	225 mm x 225 mm	10 Ltr
AUPC 252	9 x 11 Inch	225 mm x 275 mm	12 Ltr
AUPC 253	9 x 13 Inch	225 mm x 325 mm	15 Ltr
AUPC 254	12 x 12 Inch	300 mm x 300 mm	25 Ltr
AUPC 255	12 x 14 Inch	300 mm x 350 mm	26 Ltr (Discontinued)
AUPC 256	12 x 16 Inch	300 mm x 400 mm	30 Ltr
AUPC 257	12 x 18 Inch	300 mm x 450 mm	32 Ltr
AUPC 258	14 x 15 Inch	350 mm x 375 mm	39 Ltr
AUPC 259	12 x 22 Inch	300 mm x 550 mm	40 Ltr (Discontinued)
AUPC 260	14 x 17 Inch	350 mm x 425 mm	43 Ltr
AUPC 261	14 x 19 Inch	350 mm x 475 mm	45 Ltr
AUPC 262	14 x 21 Inch	350 mm x 525 mm	52 Ltr
AUPC 263	14 x 24 Inch	350 mm x 600 mm	63 Ltr