

Estd. 1966



DELIVERING TO THE WORLD

# Bacteriological Incubator Aluminium Chamber (MODEL NO : LEIN 101)



Decide with Confidence

Estd. 1966



DELIVERING TO THE WORLD

## Bacteriological Incubator Aluminium Chamber



[www.descomedicalindia.com](http://www.descomedicalindia.com)

[www.descoonline.in](http://www.descoonline.in)



# Bacteriological Incubator Aluminium Chamber (MODEL NO : LEIN 101)



## Bacteriological Incubator Aluminium Chamber Manufactured by Desco Medical

### India

#### Specification:-

Incubator has standard double wall fabrication.

Inner chamber is made out of anodized aluminium sheet and outer made out of thick mild steel sheet finished in white stoving enamel paint.

The gap is filled with special glass wool for temperature insulation.

Double wall door has double window for observation in the chamber.

Temperature ranges 5oC above ambient to 80oC + 2oC controlled by a thermostat.

L- shaped thermometer is provided for taking temperature readings.

Panel is provided with pilot lamps for mains & heaters with on-off switch and temperature Control knob.

Incubator works on 220 V AC 50 Hz single phase.

Heating is achieved from 3 sides.

Perforated removable shelves are provided.

Optional Accessories at Extra Cost :-

i) Air Circulating fan.

ii) Timer 0-24 hours.

iii) Digital temperature Indicator Cum Controller.

iv) Microprocessor PID digital temperature Indicator Cum Controller with timer

v) Digital Temperature Indicator

We are Manufacturer, Exporters and Suppliers of Bacteriological Incubator Aluminium Chamber at wholesale competitive prices from India.

Model No	Chamber Size (HXWXD) MM	Chamber Size (HXWXD) INCHES	No. of Shelves	Capacity (Approx)
LEIN 101	300 x 300 x 300 mm	12 x 12 x 12 inch	2	28 Ltrs
LEIN 102	350 x 350 x 350 mm	14 x 14 x 14 inch	2	45 Ltrs
LEIN 103	400 x 400 x 400 mm	16 x 16 x 16 inch	2	65 Ltrs
LEIN 104	450 x 450 x 450 mm	18 x 18 x 18 inch	2	95 Ltrs
LEIN 105	600 x 450 x 450 mm	24 x 18 x 18 inch	2	125 Ltrs
LEIN 106	600 x 600 x 600 mm	24 x 24 x 24 inch	2	224 Ltrs
LEIN 107	900 x 600 x 450 mm	36 x 24 x 18 inch	2	252 Ltrs
LEIN 108	900 x 600 x 600 mm	36 x 24 x 24 inch	2	336 Ltrs
LEIN 109	900 x 900 x 900 mm	36 x 36 x 36 inch	3	720 Ltrs