

Estd. 1966



DELIVERING TO THE WORLD

Blood Collection Tube Non Vacuum (MODEL NO : PPTU 601)



Blood Collection Tube Non Vacuum Manufactured by Desco Medical India



Decide with Confidence

Blood Collection Tube Non Vacuum

Estd. 1966



DELIVERING TO THE WORLD



www.descomedicalindia.com

www.descoonline.in

Specifications:-

Blood Collection Tube can provide enough and non polluted serum specimen for clinical test, while keeping the serum invariable in the long period.

Cap Colour : Red

Max. Vacuum Degree Tolerance : +10 %

Centrifuge for 10 minutes after Blood Coagulation.

It is used in Blood Collection and storage for Biochemistry, Immunology and Serology tests.

We are Manufacturer, Exporters and Suppliers of Blood Collection Tube Non Vacuum at wholesale competitive prices from India.

DELUXE SCIENTIFIC SURGICO PVT. LTD, DESCO (All Rights Reserved 2015)

Estd. 1966



DELIVERING TO THE WORLD

Blood Collection Tube Non Vacuum (MODEL NO : PPTU 601)



MODEL NO:	OUTER DIA X LENGTH	CAPACITY ML	PACK SIZE
PPTU 601	12.4 x 75 MM	2 ML	100 Pc.
PPTU 602	12.4 x 75 MM	3 ML	100 Pc.
PPTU 603	12.4 x 75 MM	4 ML	100 Pc.
PPTU 604	12.4 x 100 MM	5 ML	100 Pc.
PPTU 605	12.4 x 100 MM	6 ML	100 Pc.
PPTU 606	12.4 x 100 MM	7 ML	100 Pc.
PPTU 607	12.4 x 100 MM	8 ML	100 Pc.
PPTU 608	16 x 100 MM	9 ML	100 Pc.
PPTU 609	16 x 100 MM	10 ML	100 Pc.
PPTU 610	16 x 100 MM	11 ML	100 Pc.
PPTU 611	16 x 100 MM	12 ML	100 Pc.