

Estd. 1966



DELIVERING TO THE WORLD

Blood Collection Tube Vacuum (MODEL NO : PPTU 581)



Blood Collection Tube Vacuum Manufactured by Desco Medical India



Decide with Confidence

Blood Collection Tube Vacuum

Estd. 1966



DELIVERING TO THE WORLD



www.descomedicalindia.com

www.descoonline.in

Specifications:-

Blood Collection Tube can provide enough and non polluted serum specimen for clinical test, while keeping the serum invariable in the long period.

Cap Colour : Red

Max. Vacuum Degree Tolerance : +10 %

Centrifuge for 10 minutes after Blood Coagulation.

It is used in Blood Collection and storage for Biochemistry, Immunology and Serology tests.

We are Manufacturer, Exporters and Suppliers of Blood Collection Tube Vacuum at wholesale competitive prices from India.

Estd. 1966



DELIVERING TO THE WORLD

Blood Collection Tube Vacuum (MODEL NO : PPTU 581)



MODEL NO:	OUTER DIA X LENGTH	CAPACITY ML	PACK SIZE
PPTU 581	12.4 x 75 MM	2 ML	100 Pc.
PPTU 582	12.4 x 75 MM	3 ML	100 Pc.
PPTU 583	12.4 x 75 MM	4 ML	100 Pc.
PPTU 584	12.4 x 100 MM	5 ML	100 Pc.
PPTU 585	12.4 x 100 MM	6 ML	100 Pc.
PPTU 586	12.4 x 100 MM	7 ML	100 Pc.
PPTU 587	12.4 x 100 MM	8 ML	100 Pc.
PPTU 588	16 x 100 MM	9 ML	100 Pc.
PPTU 589	16 x 100 MM	10 ML	100 Pc.
PPTU 590	16 x 100 MM	11 ML	100 Pc.
PPTU 591	16 x 100 MM	12 ML	100 Pc.