

Estd. 1966



DELIVERING TO THE WORLD

Cautery Machine 300 Watt Analog(MODEL NO : IECM 101)



Cautery Machine 300 Watt Analog Supplied by Desco Medical India



Decide with Confidence

Cautery Machine 300 Watt Analog

Estd. 1966



DELIVERING TO THE WORLD



www.descomedicalindia.com

www.descoonline.in

Estd. 1966



Cautery Machine 300 Watt Analog(MODEL NO : IECM 101)



Specification:-

Cautery Machine is also known as Electrocautery Machine

Electrocautery is used in minor surgical procedures in dermatology, ophthalmology, otolaryngology, plastic surgery, and urology.

In Electrocautery procedure, heated electrode is applied to living tissue to achieve hemostasis or varying degrees of tissue destruction.

300Watt analog model with basic cut, coagulation and bipolar modes with capability of doing surgeries.

Power : 300 W Analog Non Programmable

Capability : basic surgeries Unbelievable economy for basic functionality

Max Power : 300W

Programs : NA

Monopolar Modes : Cut, Coagulation

Bipolar Modes : Cut, Coag

Tissue Response Circuit : NA

Patient Plate Monitoring System : NA

Dimension (W x H x D) & Weight : 3.5 " X 13 X" & 6 Kg.

Standard accessories with the unit :-

1 Surgical Monopolar active handle valley lab type. 1nos

2 Needle set (electrodes set for above probe) 1nos

3 S.S.patient plate (adult) 1nos

4 Patient plate cable cord. 1nos

5 Hand switch pencil disposable.(1time use) 2nos

6 Double paddle footswitch. 1nos

7 Bipolar forceps regular bayonet. 1nos

8 Bipolar forceps cable cord. 1nos

We are Exporters and Suppliers of Cautery Machine 300 Watt Analog at wholesale competitive prices from India.

MODEL NO

IECM 101

Cautery Macine 300 Watt Analog