

Estd. 1966



DELIVERING TO THE WORLD

Cautery Machine 300 Watt Digital(MODEL NO : IECM 105)



Cautery Machine 300 Watt Digital Supplied by Desco Medical India



Decide with Confidence

Cautery Machine 300 Watt Digital

Estd. 1966



DELIVERING TO THE WORLD



www.descomedicalindia.com

www.descoonline.in

Estd. 1966



DELIVERING TO THE WORLD

Cautery Machine 300 Watt Digital(MODEL NO : IECM 105)



Specification:-

Cautery Machine is also known as Electrocautery Machine and Electro Surgical unit. Electrocautery is used in minor surgical procedures in dermatology, ophthalmology, otolaryngology, plastic surgery, and urology.

In Electrocautery procedure, heated electrode is applied to living tissue to achieve hemostasis or varying degrees of tissue destruction.

Cautery Machine 300 Watt Digital is micro controller based model with 6 programs

Power Digital : 300 W

Programs : 6

Capability : All surgeries but till 300 W power

Ideal where 400 W not required

Max Power : 300W

Display : 7 Segment Digital

Programs : 6

Monopolar Cut Modes : Pure Cut, Blend1, Blend2, Endo Cut

Monopolar Coagulation Modes : Spray, Dessicate, Fulgurate

Bipolar Modes : Micro, Macro

Tissue Response Circuit : Present

Patient Plate Monitoring System : Present

Dimension (WxHxD) & Weight : 4.5" x 10" x 11" & 3.8Kg.

Standard accessories with the unit :-

1 Surgical Monopolar active handle valley lab type. 1nos

2 Needle set (electrodes set for above probe) 1nos

3 S.S.patient plate (adult) 1nos

4 Patient plate cable cord. 1nos

5 Hand switch pencil disposable.(1time use) 2nos

6 Double paddle footswitch. 1nos

7 Bipolar forceps regular bayonet. 1Nos

8 Bipolar forceps cable cord. 1nos

9 Single paddle foot switch 1nos

We are Exporters and Suppliers of Cautery Machine 300 Watt Digital at wholesale competitive prices from India.

MODEL NO

IECM 105

Cautery Macine 300 Watt Digital