

Estd. 1966



DELIVERING TO THE WORLD

Cautery Machine 400 Watt Analog(MODEL NO : IECM 110)



Cautery Machine 400 Watt Analog Supplied by Desco Medical India



Decide with Confidence

Cautery Machine 400 Watt Analog

Estd. 1966



DELIVERING TO THE WORLD



www.descomedicalindia.com

www.descoonline.in

Estd. 1966



Cautery Machine 400 Watt Analog(MODEL NO : IECM 110)



Specification:-

Cautery Machine is also known as Electrocautery Machine

Electrocautery is used in minor surgical procedures in dermatology, ophthalmology, otolaryngology, plastic surgery, and urology.

In Electrocautery procedure, heated electrode is applied to living tissue to achieve hemostasis or varying degrees of tissue destruction.

Cautery Machine is a 400W analog model with number of cut, blend and coagulation modes to do under water surgery.

Power : 400 W Analog Non Programmable

Capability : All surgeries

Robust quality : Great price

Max Power : 400W

Programs : NA

Monopolar Cut Modes : Pure Cut, Blend1, Blend2

Monopolar Coagulation Modes : Spray, Dessicate, Fulgurate

Bipolar Modes : Coag

Tissue Response Circuit : NA

Patient Plate Monitoring System : Present

Dimension (WxHxD) & Weight : 4.5" x 10" x 11" & 3.8 Kg.

Standard accessories with the unit :-

1 Surgical Monopolar active handle valley lab type. 1nos

2 Needle set (electrodes set for above probe) 1nos

3 S.S.patient plate (adult) 1nos

4 Patient plate cable cord. 1nos

5 Hand switch pencil disposable.(1time use) 2nos

6 Double paddle footswitch. 1nos

7 Bipolar forceps regular bayonet. 1nos

8 Bipolar forceps cable cord. 1nos

9 Single paddle foot switch 1nos

We are Exporters and Suppliers of Cautery Macine 400 Watt Analog at wholesale competitive prices from India.

MODEL NO

IECM 110

Cautery Macine 400 Watt Analog