

Estd. 1966



DELIVERING TO THE WORLD

Clot Activator Tube Non Vacuum (MODEL NO : PPTU 521)



Clot Activator Tube Non Vacuum Manufactured by Desco Medical India



Decide with Confidence

Clot Activator Tube Non Vacuum



DELIVERING TO THE WORLD



www.descomedicalindia.com

www.descoonline.in

Specifications:-

Clot Activator Tube is Coated with micronised Silica particles which causes blood to clot rapidly. Clot Activator Tube provides enough and non polluted serum specimen for clinical test.

Cap Colour : Red

Max. Vacuum Degree Tolerance : +10 %

Invert the Tube for 6-8 times gently and centrifuge for 10 minutes after 15-30 minutes coagulation.

It is used in Blood Collection and storage for Biochemistry, Immunology and Serology tests.

We are Manufacturer, Exporters and Suppliers of Clot Activator Tube Non Vacuum at wholesale competitive prices from India.

Estd. 1966



Clot Activator Tube Non Vacuum (MODEL NO : PPTU 521)



DELIVERING TO THE WORLD

MODEL NO:	OUTER DIA X LENGTH	CAPACITY ML	PACK SIZE
PPTU 521	12.4 x 75 MM	2 ML	100 Pc.
PPTU 522	12.4 x 75 MM	3 ML	100 Pc.
PPTU 523	12.4 x 75 MM	4 ML	100 Pc.
PPTU 524	12.4 x 100 MM	5 ML	100 Pc.
PPTU 525	12.4 x 100 MM	6 ML	100 Pc.
PPTU 526	12.4 x 100 MM	7 ML	100 Pc.
PPTU 527	12.4 x 100 MM	8 ML	100 Pc.
PPTU 528	16 x 100 MM	9 ML	100 Pc.
PPTU 529	16 x 100 MM	10 ML	100 Pc.
PPTU 530	16 x 100 MM	11 ML	100 Pc.
PPTU 531	16 x 100 MM	12 ML	100 Pc.