

Continuous Passive Motion Unit (L/E) Supplied by Desco Medical India


CE


D&B
Decide with Confidence

Estd. 1966


DESCO[®]
DELIVERING TO THE WORLD

Continuous Passive Motion Unit (L/E)


www.descomedicalindia.com



www.descoonline.in

Specification:-

Knee CPM, a micro compute controlled unit, represents the latest breakthrough in CPM technology.

Solid state circuitry control panel with highly visual LED Digital readouts allows flexion and extension of knee joint through a prefixed range of motion and time.

Proper fit of calf and thigh lengths are easily achieved.

The necessary limb support provide comfortable surface for prolonged contact & are washable in cold water.

The light weight units is constructed of aluminium sheet & is place on a movable trolley.



Continuous Passive Motion Unit (L/E) (MODEL NO : PTCU 101)



Continuous Passive Motion Unit (L/E) Supplied by Desco Medical India

Specification:-

Technical :

- * Therapy Modes : Six, Mobilization Modes for continuous, Intermittent and Progressive passive movements.
- * Rom Setting : Programmable from 00 to 110 degree.
- * Pause Time : Programmable Extension / Flexion Paused from 00 to 99 seconds.
- * Treatment Time : Programmable from 01 to 99 min.
- * Inching Switch : Starting Angle can be set according to the need of patient.
- * Memory / Programme : 30 Treatment Protocols can be programmed & stored in memory.
- * Displays : Digital display of ROM setting Pause Time, Treatment Time & Graphic display of Mobilization modes.
- * Patient Safety : Patient Stop / interrupt treatment with hand held switch if needed.
- * Open Programs : 20 Programs are free to store your own patterns.
- * Trolley : Provides Transportation & Storage.
- * Power Supply : 200-240V AC at 50-60 Hz.

We are Exporters and Suppliers of Continuous Passive Motion Unit (L/E) at wholesale competitive prices from India.

MODEL NO:

PTCU 101 Continuous Passive Motion Unit (L/E)