

Estd. 1966



DELIVERING TO THE WORLD

Corrugated Breathing Tube (MODEL NO : APAA 131)



Decide with Confidence

Corrugated Breathing Tubes

Estd. 1966



DELIVERING TO THE WORLD



www.descomedicalindia.com

www.anesthesiaproductsindia.com

www.descoonline.in

Corrugated Breathing Tube Manufactured by Desco India

Specification:-

Corrugated Breathing Tube is made of high grade antistatic latex rubber.

The inner dia size is 15 mm & 22 mm as per international standard.

Corrugated Breathing Tube is a flexible tube for breathing through

It can be used as part of a scuba set or other breathing apparatus or a medical oxygen apparatus or anaesthetic apparatus

These breathing tubes are also distinguished from the medium-pressure hoses which are often found as parts of modern breathing apparatus

They are wide, and usually corrugated to let the user's head move about without the tube pinching at kinks.

DELUXE SCIENTIFIC SURGICO PVT. LTD, DESCO (All Rights Reserved 2015)

Estd. 1966



DELIVERING TO THE WORLD

Corrugated Breathing Tube (MODEL NO : APAA 131)



Each end has a screw connection.

They may contain a one-way valve to keep the air or gas flowing the right way.

We are Manufacturer, Exporters and Suppliers of Corrugated Breathing Tube at wholesale competitive prices from India.

Model No:

APAA 131 12 Inch

APAA 132 16 Inch

APAA 133 42 Inch