

# Dispovan Needles (MODEL NO : MDSN 101)



## Dispovan Needles Manufactured by Desco Medical India



### Specification:-

Made from imported cold rolled stainless steel strips.

Tips of the needles are ground on high precision CNC machines.

Siliconised, minimizing drag during penetration and withdrawal.

Needle point is checked by vision camera system.

#### COMPONENTS:

Needle Hubs are made of non-toxic, medical-grade polypropylene compatible with any injectable medication.

Cannulas Burr Free 3 Facet Lancet Point are made from SUS stainless steel.

#### Packing :

The needles are packed in tamper-evident, blister packs having easily peelable, fibre-free, imported, medical-grade paper to ensure reliability of sealing and product sterility.

We are Exporters and Suppliers of Dispovan Needles at wholesale competitive prices from India.



# Dispovan Needles (MODEL NO : MDSN 101)



<b>Model No</b>	<b>Gauge Size</b>	<b>Outer Dia in mm</b>	<b>Colour Code</b>	<b>Length in mm</b>
MDSN 101	16 G	1.60	Natural	25 mm
MDSN 102	16 G	1.60	Natural	38 mm
MDSN 103	18 G	1.20	Pink	18 mm
MDSN 104	18 G	1.20	Pink	25 mm
MDSN 105	18 G	1.20	Pink	38 mm
MDSN 106	19 G	1.10	Cream	32 mm
MDSN 107	19 G	1.10	Cream	38 mm
MDSN 108	20 G	0.90	Yellow	25 mm
MDSN 109	20 G	0.90	Yellow	32 mm
MDSN 110	20 G	0.90	Yellow	38 mm
MDSN 111	21 G	0.80	Green	25 mm
MDSN 112	21 G	0.80	Green	32 mm
MDSN 113	21 G	0.80	Green	38 mm
MDSN 114	22 G	0.70	Black	18 mm
MDSN 115	22 G	0.70	Black	25 mm
MDSN 116	22 G	0.70	Black	32 mm
MDSN 117	22 G	0.70	Black	38 mm
MDSN 118	23 G	0.60	Blue	25 mm
MDSN 119	23 G	0.60	Blue	32 mm
MDSN 120	23 G	0.60	Blue	38 mm
MDSN 121	24 G	0.55	Purple	25 mm
MDSN 122	24 G	0.55	Purple	32 mm
MDSN 123	24 G	0.55	Purple	38 mm
MDSN 124	25 G	0.50	Orange	16 mm
MDSN 125	25 G	0.50	Orange	25 mm
MDSN 126	25 G	0.50	Orange	38 mm
MDSN 127	26 G	0.45	Brown	4 mm
MDSN 128	26 G	0.45	Brown	13 mm
MDSN 129	26 G	0.45	Brown	38 mm
MDSN 130	27 G	0.40	Grey	13 mm
MDSN 131	27 G	0.40	Grey	18 mm