

## Dispovan Syringe with Needle Manufactured by Desco Medical India

**CE**  **Decide with Confidence**

**Dispovan Syringe with Needle**

**DESCO**<sup>®</sup>  
Estd. 1966  
DELIVERING TO THE WORLD



[www.descomedicalindia.com](http://www.descomedicalindia.com) [www.descoonline.in](http://www.descoonline.in)

### Specification:-

Luer mount Single use syringes with Needle.

Barrels are made of nontoxic, medical-grade polypropylene compatible with any medication.

The tip of the barrel is luer mount type having a 6% luer taper.

Plungers are made of nontoxic, medical-grade polypropylene compatible with any medication.

Gaskets are made of natural rubber which is chemically inert and compatible with medication for short term contact drug delivery application.

We are Exporters and Suppliers of Dispovan Syringe with Needle at wholesale competitive prices from India.



# Dispovan Syringe with Needle (MODEL NO : MDSN 451)



<b>Model No</b>	<b>Capacity</b>	<b>Gauge</b>	<b>Length in mm</b>	<b>Dia in mm</b>	<b>Type of Nozzle</b>
MDSN 451	2 ml	21G	25 mm	0.80	Luer Mount with Central Nozzle
MDSN 452	2 ml	21G	38 mm	0.80	Luer Mount with Central Nozzle
MDSN 453	2 ml	22G	25 mm	0.70	Luer Mount with Central Nozzle
MDSN 454	2 ml	22G	32 mm	0.70	Luer Mount with Central Nozzle
MDSN 455	2 ml	22G	38 mm	0.70	Luer Mount with Central Nozzle
MDSN 456	2 ml	23G	25 mm	0.60	Luer Mount with Central Nozzle
MDSN 457	2 ml	23G	38 mm	0.60	Luer Mount with Central Nozzle
MDSN 458	2 ml	24G	25 mm	0.55	Luer Mount with Central Nozzle
MDSN 459	3 ml	21G	25 mm	0.80	Luer Mount with Central Nozzle
MDSN 460	3 ml	21G	38 mm	0.80	Luer Mount with Central Nozzle
MDSN 461	3 ml	22G	25 mm	0.70	Luer Mount with Central Nozzle
MDSN 462	3 ml	22G	32 mm	0.70	Luer Mount with Central Nozzle
MDSN 463	3 ml	22G	38 mm	0.70	Luer Mount with Central Nozzle
MDSN 464	3 ml	23G	25 mm	0.60	Luer Mount with Central Nozzle
MDSN 465	3 ml	23G	38 mm	0.60	Luer Mount with Central Nozzle
MDSN 466	3 ml	24G	25 mm	0.55	Luer Mount with Central Nozzle
MDSN 467	5 ml	21G	25 mm	0.80	Luer Mount with Central Nozzle
MDSN 468	5 ml	21G	38 mm	0.80	Luer Mount with Central Nozzle
MDSN 469	5 ml	22G	25 mm	0.70	Luer Mount with Central Nozzle
MDSN 470	5 ml	22G	32 mm	0.70	Luer Mount with Central Nozzle
MDSN 471	5 ml	22G	38 mm	0.70	Luer Mount with Central Nozzle
MDSN 472	5 ml	23G	25 mm	0.60	Luer Mount with Central Nozzle
MDSN 473	5 ml	23G	38 mm	0.60	Luer Mount with Central Nozzle
MDSN 474	5 ml	24G	25 mm	0.55	Luer Mount with Central Nozzle
MDSN 475	10 ml	20G	38 mm	0.90	Luer Mount with Eccentric Nozzle
MDSN 476	10 ml	21G	38 mm	0.80	Luer Mount with Eccentric Nozzle
MDSN 477	10 ml	22G	32 mm	0.70	Luer Mount with Eccentric Nozzle
MDSN 478	10 ml	22G	38 mm	0.70	Luer Mount with Eccentric Nozzle
MDSN 479	20 ml	20G	38 mm	0.90	Luer Mount with Eccentric Nozzle
MDSN 480	20 ml	21G	38 mm	0.80	Luer Mount with Eccentric Nozzle
MDSN 481	20 ml	22G	32 mm	0.70	Luer Mount with Eccentric Nozzle
MDSN 482	20 ml	22G	38 mm	0.70	Luer Mount with Eccentric Nozzle