

Estd. 1966



DELIVERING TO THE WORLD

EDTA K3 Tube Vacuum (MODEL NO : PPTU 741)



EDTA K3 Tube Vacuum Manufactured by Desco Medical India



Decide with Confidence

EDTA K3 Tube Vacuum



Estd. 1966
DELIVERING TO THE WORLD



www.descomedicalindia.com www.descoonline.in

Specifications:-

EDTA K3 Tubes can effectively chelat the calcium ion of blood which can prevent the blood \coagulation.

EDTA K3 Tube does not affect the count and size of the leukocyte and keep erythrocyte invariable.

Cap Colour : Lavender or Purple

Max. Vacuum Degree Tolerance : +10 %

Invert the Tube for 6-8 times immediately after drawing blood.

EDTA K3 Tube is used in Clinical hematology test & blood banking.

We are Manufacturer, Exporters and Suppliers of EDTA K3 Tube Vacuum at wholesale competitive prices from India.

Estd. 1966



DELIVERING TO THE WORLD

EDTA K3 Tube Vacuum (MODEL NO : PPTU 741)



MODEL NO:	OUTER DIA X LENGTH	CAPACITY ML	PACK SIZE
PPTU 741	12.4 x 75 MM	2 ML	100 Pc.
PPTU 742	12.4 x 75 MM	3 ML	100 Pc.
PPTU 743	12.4 x 75 MM	4 ML	100 Pc.
PPTU 744	12.4 x 100 MM	5 ML	100 Pc.
PPTU 745	12.4 x 100 MM	6 ML	100 Pc.
PPTU 746	12.4 x 100 MM	7 ML	100 Pc.
PPTU 747	12.4 x 100 MM	8 ML	100 Pc.
PPTU 748	16 x 100 MM	9 ML	100 Pc.
PPTU 749	16 x 100 MM	10 ML	100 Pc.
PPTU 750	16 x 100 MM	11 ML	100 Pc.
PPTU 751	16 x 100 MM	12 ML	100 Pc.