

Estd. 1996



# ELECTRIC C-ARM TABLE WITH OVERRIDE PANEL (MODEL NO : OTCE 120)



DELIVERING TO THE WORLD

## Electric C-Arm Table With Override Panel Manufactured by Desco Medical India

  
**CE**

  
BSCIC  
ISO 9001 REGISTERED

  
D&B  
Decide with Confidence

  
Estd. 1966  
DESCO  
DELIVERING TO THE WORLD

### ELECTRIC C-ARM TABLE WITH OVERRIDE PANEL



[www.descomedicalindia.com](http://www.descomedicalindia.com)[www.operatingtablesindia.com](http://www.operatingtablesindia.com)

### Specification:-

- Table Top Size (LxW) : 2150 x 535 mm
- Minimum Height (Without) : 680 mm
- Maximum Height (Without) : 1000 mm
- Trendelenburg /Reverse Trend : 45 degree/30 degree
- Lateral Tilt (Left/Right) : 20 degree
- Back Section adjustment : + 70 degree/-35 degree
- Leg Section adjustment : 0 degree / -90 degree & Split
- Longitudinal Slide : 300 mm
- Battery Backup : 24 hours
- Power Supply : 220VAC -50Hz
- Load Carrying Capacity : 350 Kg



# ELECTRIC C-ARM TABLE WITH OVERRIDE PANEL (MODEL NO : OTCE 120)



DELIVERING TO THE WORLD

- Six section eccentrically positioned radio-translucent table top without cross members facilitates C-arm imaging.
  - Operating Positions i.e. Height, Trendelenburg, Reverse Trendelenburg, Lateral Tilt, Back Section
  - Longitudinal slide are quickly & precisely controlled by handheld remote as well as a Backup remote with symbolic figures mounted on the table.
  - Imported Electro-Hydraulic System (From Netherlands) comprising of Hydraulic Motor pump, Solenoid Valves & a manual Override System to assure smooth & quiet movements.
  - Battery Backup: The table has built-in Battery Backup System with enough power for several hours of operations in case of power failure.
  - Manual Override System : The table is equipped with a foot pedal operated Manual Override System with a position selector for height, Trendelenburg, Lateral Tilt & Back Section adjustment. This can be used in case of Handset failure, Power failure, Battery failure. The manual override system works independently and is not linked with any electric function.
  - Built-in Kidney Elevator (C-Arm Compatible)- manually operated.
- State of art Maintenance Free Self Lubricating Lifting system with hi-precision bearings with large contact surface assures maximum stability.
- Non-Corrosive 55304 Grade Stainless Steel Base & Cylinders covers, side rails, clamps & Standard Accessories for easy clean in g and hygiene.
- Bariatric Surgery Ready: This table can be used for Bariatric surgeries by using width extenders. Table capable of achieving extreme Trendelenburg of 45 Degrees. Load carrying capacity of 350 Kg.

Standard Accessories :-

Shaped anesthetic Screen : 1 Pc

Padded shoulder support : 1 Pair

Padded lateral support : 1 Pair

Radio translucent armboard with pad :1 Pair

Geopal Knee crutch with straps : 1 Pair

PU Foam Matters : 1 Set

Optional Accessories :-

Side Extensions for Bariatric Surgery with body Strap

Foot Rest : 1 Pair

We are Manufacturer, Exporters and Suppliers of Electric C-Arm Table With Override Panel at wholesale competitive prices from India.

## **MODEL NO:**

**OTCE 120**      Electric C-Arm Table With Override Panel