

Estd. 1966



# Electric Luxury Obstetric Table (MODEL NO : DTEL 105)



DELIVERING TO THE WORLD

## Electric Luxury Obstetric Table Manufactured by Desco Medical India



### Specification:-

This Table is Suitable for Examination Diagnosis & Operation for Gynaecology, Urological Surgery.

### Operating Positions

Three Sectional upholstered tabletop mounted on strong S.S. tubular frame ensures maximum comfort to the patient.

- Maintenance free lifting system with heavy duty ball bearings provide silent & jerk free movements.
- Built-in battery backup for emergency use only.
- Fully Electric Table with Operating Positions.

### Height Adjustment

- Trendelenburg & Reverse Trendelenburg.
- Back Section Lift.
- SS 304 Grade steel Used.
- Leg section is easily detachable and floor locking is manual.
- Table is suitable for examination diagnosis and operation for Gynecology, Urological Surgery.
- Actuators used Ti Motion (Taiwan).
- **Dual Access control**-Hand held Remote & Foot Switch.
- Specially design base with 'U' Cut providing extra leg room while examination.

Estd. 1966



# Electric Luxury Obstetric Table (MODEL NO : DTEL 105)



DELIVERING TO THE WORLD

- \* Table Size (LXW) : 191 x 62 cm.
- \* Min. Height – 78 cm.
- \* Max. Height – 102cm.
- \* Trendelenburg / Rev. Trendelenburg – 25 Degree.
- \* Back Section - 70 Degree.

## **Standard Accessories**

- \* 'L' Shaped Anaesthetic Screen : 1 Pc.
- \* Padded Shoulder Support : 1 Pair.
- \* Radio-translucent Armboard with pad : 2 Pcs.
- \* Goepel Knee Crutch with straps : 1 Pair.
- \* PU Foam Mattress : 1 Set.
- \* Water proof table top cover : 1 Pc.
- \* S.S. Douche Tray with drain pipe : 1 Pc.

We are Manufacturer, Exporters and Suppliers of Electric Luxury Obstetric Table at wholesale competitive prices from India.

## **MODEL NO:**

**DTEL 105**

Electric Luxury Obstetric Table