

Estd. 1966



DELIVERING TO THE WORLD

# Fluoride Tube Vacuum (MODEL NO : PPTU 781)



## Fluoride Tube Vacuum Manufactured by Desco Medical India



Decide with Confidence

### Fluoride Tube Vacuum

Estd. 1966



DELIVERING TO THE WORLD



[www.descomedicalindia.com](http://www.descomedicalindia.com)

[www.descoonline.in](http://www.descoonline.in)

### Specifications:-

Fluoride Tube has additive EDTA potassium salt

EDTA potassium salt can chelate the calcium ion of blood prevents the blood coagulation.

Sodium Fluoride control the glycolysis.

Cap Colour : Grey

Max. Vacuum Degree Tolerance : +10 %

Invert the Tube for 6-8 times gently.

Fluoride Tube is used for the analysis of blood glucose & lactate.

We are Manufacturer, Exporters and Suppliers of Fluoride Tube Vacuum at wholesale competitive prices from India.

Estd. 1966



DELIVERING TO THE WORLD

# Fluoride Tube Vacuum (MODEL NO : PPTU 781)



<b>MODEL NO:</b>	<b>OUTER DIA X LENGTH</b>	<b>CAPACITY ML</b>	<b>PACK SIZE</b>
PPTU 781	12.4 x 75 MM	2 ML	100 Pc.
PPTU 782	12.4 x 75 MM	3 ML	100 Pc.
PPTU 783	12.4 x 75 MM	4 ML	100 Pc.
PPTU 784	12.4 x 100 MM	5 ML	100 Pc.
PPTU 785	12.4 x 100 MM	6 ML	100 Pc.
PPTU 786	12.4 x 100 MM	7 ML	100 Pc.
PPTU 787	12.4 x 100 MM	8 ML	100 Pc.
PPTU 788	16 x 100 MM	9 ML	100 Pc.
PPTU 789	16 x 100 MM	10 ML	100 Pc.
PPTU 790	16 x 100 MM	11 ML	100 Pc.
PPTU 791	16 x 100 MM	12 ML	100 Pc.