

Estd. 1966



DELIVERING TO THE WORLD

Gel & Clot Activator Tube Non Vacuum (MODEL NO : PPTU 561)



Gel & Clot Activator Tube Non Vacuum Manufactured by Desco Medical India



Specifications:-

Gel & Clot Activator Tube is Coated with micronised Silica particles which causes blood to clot rapidly and the barrier gel can effectively separate the serum from fibrin and cells, while preventing substance exchange between blood cell and serum.

Owing to its clot activator, it is better for keeping the original character of the blood and solving the problem of hemolysis.

Cap Colour : Golden Yellow

Max. Vacuum Degree Tolerance : +10 %

Invert the Tube for 6-8 times gently and centrifuge for 10 minutes after 3000-3500rpm after 15-30 minutes coagulation.

It is used in Blood Collection and storage for Biochemistry, Immunology and Serology tests.

We are Manufacturer, Exporters and Suppliers of Gel & Clot Activator Tube Non Vacuum at wholesale competitive prices from India.

DELUXE SCIENTIFIC SURGICO PVT. LTD, DESCO (All Rights Reserved 2015)

Estd. 1966



DELIVERING TO THE WORLD

Gel & Clot Activator Tube Non Vacuum (MODEL NO : PPTU 561)



MODEL NO:	OUTER DIA X LENGTH	CAPACITY ML	PACK SIZE
PPTU 561	12.4 x 75 MM	2 ML	100 Pc.
PPTU 562	12.4 x 75 MM	3 ML	100 Pc.
PPTU 563	12.4 x 75 MM	4 ML	100 Pc.
PPTU 564	12.4 x 100 MM	5 ML	100 Pc.
PPTU 565	12.4 x 100 MM	6 ML	100 Pc.
PPTU 566	12.4 x 100 MM	7 ML	100 Pc.
PPTU 567	12.4 x 100 MM	8 ML	100 Pc.
PPTU 568	16 x 100 MM	9 ML	100 Pc.
PPTU 569	16 x 100 MM	10 ML	100 Pc.
PPTU 570	16 x 100 MM	11 ML	100 Pc.
PPTU 571	16 x 100 MM	12 ML	100 Pc.