



Hot Air Universal Oven Aluminium Chamber (MODEL NO : LEOV 101)



Decide with Confidence

Hot Air Universal Oven Aluminium Chamber

Estd. 1966



DELIVERING TO THE WORLD



www.descomedicalindia.com

www.descoonline.in

Hot Air Universal Oven Aluminium Chamber Manufactured by Desco Medical India

Specification:-

Hot Air Oven Universal are also known as Hot Air Oven Memmert type

Inner chamber is made of anodized aluminum sheet and outer is made of thick mild steel sheet finished in white stoving enamel paint.

Oven has double wall fabrication. The gap between the walls is filled with special grade glass wool for heat insulation.

Inner chamber is ribbed for placing two or three removable stainless steel shelves at convenient levels.

Ovens have three side heating elements made of nichrome wire to attain quick and uniform heating in range of 50oC to 250oC + 1oC controlled by capillary type thermostat.

L- shaped thermometer is provided for taking temperature readings.

Control Panel is provided with selector switch of high or low rates of power, temperature Control knob and indicators for mains & thermostat

Incubator works on 220 V AC 50 Hz single phase.

Optional Accessories at Extra Cost :-

- i) Air Circulating fan.
- ii) Timer 0-24 hours.
- iii) Digital temperature Indicator Cum Controller.
- iv) Microprocessor PID digital temperature Indicator Cum Controller with timer
- v) Digital Temperature Indicator

We are Manufacturer, Exporters and Suppliers of Hot Air Universal Oven Aluminium Chamber at wholesale competitive prices from India.

Model No	Chamber Size (HXWxD) MM	Chamber Size (HXWxD) INCHES	No. of Shelves	Capacity (Approx)
LEOV 101	300 x 300 x 300 mm	12 x 12 x 12 inch	2	28 Ltrs
LEOV 102	350 x 350 x 350 mm	14 x 14 x 14 inch	2	45 Ltrs
LEOV 103	400 x 400 x 400 mm	16 x 16 x 16 inch	2	65 Ltrs
LEOV 104	450 x 450 x 450 mm	18 x 18 x 18 inch	2	95 Ltrs
LEOV 105	600 x 450 x 450 mm	24 x 18 x 18 inch	2	125 Ltrs
LEOV 106	600 x 600 x 600 mm	24 x 24 x 24 inch	2	224 Ltrs
LEOV 107	900 x 600 x 450 mm	36 x 24 x 18 inch	2	252 Ltrs
LEOV 108	900 x 600 x 600 mm	36 x 24 x 24 inch	2	336 Ltrs
LEOV 109	900 x 900 x 900 mm	36 x 36 x 36 inch	3	720 Ltrs