

Estd. 1966



DELIVERING TO THE WORLD

Hot Air Universal Oven Stainless Steel Chamber (MODEL NO : LEOV 151)



Decide with Confidence

Estd. 1966



DELIVERING TO THE WORLD

Hot Air Universal Oven Stainless Steel Chamber



www.descomedicalindia.com

www.descoonline.in



Hot Air Universal Oven Stainless Steel Chamber (MODEL NO : LEOV 151)



Hot Air Universal Oven Stainless Steel Chamber Manufactured by Desco Medical India

Specification:-

Hot Air Oven Universal are also known as Hot Air Oven Memmert type

Inner chamber is made of stainless sheet and outer is made of thick mild steel sheet finished in white stoving enamel paint.

Oven has double wall fabrication. The gap between the walls is filled with special grade glass wool for heat insulation.

Inner chamber is ribbed for placing two or three removable stainless steel shelves at convenient levels.

Ovens have three side heating elements made of nichrome wire to attain quick and uniform heating in range of 50oC to 250oC + 1oC controlled by capillary type thermostat.

L- shaped thermometer is provided for taking temperature readings.

Control Panel is provided with selector switch of high or low rates of power, temperature

Control knob and indicators for mains & thermostat

Incubator works on 220 V AC 50 Hz single phase.

Optional Accessories at Extra Cost :-

i) Air Circulating fan.

ii) Timer 0-24 hours.

iii) Digital temperature Indicator Cum Controller.

iv) Microprocessor PID digital temperature Indicator Cum Controller with timer

v) Digital Temperature Indicator

We are Manufacturer, Exporters and Suppliers of Hot Air Universal Oven Stainless Steel Chamber at wholesale competitive prices from India.

Model No	Chamber Size (HXWxD) MM	Chamber Size (HXWxD) INCHES	No. of Shelves	Capacity (Approx)
LEOV 151	300 x 300 x 300 mm	12 x 12 x 12 inch	2	28 Ltrs
LEOV 152	350 x 350 x 350 mm	14 x 14 x 14 inch	2	45 Ltrs
LEOV 153	400 x 400 x 400 mm	16 x 16 x 16 inch	2	65 Ltrs
LEOV 154	450 x 450 x 450 mm	18 x 18 x 18 inch	2	95 Ltrs
LEOV 155	600 x 450 x 450 mm	24 x 18 x 18 inch	2	125 Ltrs
LEOV 156	600 x 600 x 600 mm	24 x 24 x 24 inch	2	224 Ltrs
LEOV 157	900 x 600 x 450 mm	36 x 24 x 18 inch	2	252 Ltrs
LEOV 158	900 x 600 x 600 mm	36 x 24 x 24 inch	2	336 Ltrs
LEOV 159	900 x 900 x 900 mm	36 x 36 x 36 inch	3	720 Ltrs