

Estd. 1966



DELIVERING TO THE WORLD

Hydraulic OT Table with Side End Control (MODEL NO : OTHY 105)



Hydraulic OT Table with Side End Control Manufactured by Desco Medical India



Decide with Confidence

Hydraulic OT Table with Side End Control

Estd. 1966



DELIVERING TO THE WORLD



INCLUDED Standard Accessories



www.descomedicalindia.com

www.descoonline.in

www.operatingtablesindia.com

Specification:-

- * Stainless Steel Fitting.
- * Detachable head, leg & Pelvic section.
- * Base & Cover made of attractive, easy to clean Stainless Steel.
- * Sophisticated mechanics provide smooth, step less articulation of the table Top for precise patient positioning.
- * Center Pillar / In Built Kidney Bridge.

Estd. 1966



DELIVERING TO THE WORLD

Hydraulic OT Table with Side End Control (MODEL NO : OTHY 105)



Technical Specs:-

Length of Table	1900 mm
Width of Table Top	520 mm
Minimum Height (without Mattress)	Min. 750 mm
Maximum Height (without Mattress)	Max. 1000 mm
Trendelenburg	25°
Reverse Trendelenburg	25°
Lateral Tilt	±20°
Back Section	±60°
Head Section	±60°
Leg Section	0°-90°
Floor Locking System	Manual

We are Manufacturer, Exporters and Suppliers of Hydraulic OT Table with Side End Control at wholesale competitive prices from India.

MODEL NO:

OTHY 105

Hydraulic OT Table with Side End Control