

Infant Radiant Warmer Diamond with Phototherapy manufactured by Desco Medical India


CE
Decide with Confidence


D&B
Decide with Confidence

Estd. 1966


DESCO[®]
DELIVERING TO THE WORLD

Infant Radiant Warmer Diamond With Phototherapy

CFL Phototherapy





LED Phototherapy



www.descomedicalindia.com

www.babycareequipments.com

www.descoonline.in

Specification:-

Infant Radiant Warmer is digital micro processor - based controller comes with a full spectrum of features.

Control Panel : The Infant Radiant Warmer comes complete with automatic manual /baby (Servo) / Air modes.

In addition, the Skin/Air temperature display is larger than set temperature display.

Skin temperature display can be converted form OC to OF by the push of a button.

Two function control modes provide more flexibility while optimizing the thermal environment.

A big LCD is provided for more information and visualizations of the heater power in manual and servo modes.

A full complement of alarms provide safety and convenience, with both audible and visual indicators in the event of :



Infant Radiant Warmer Diamond With Phototherapy (MODEL NO : BCRW 121)



- High Temp./ Low Temp.
- System Fail
- Probe Failure
- Power Failure
- Timer
- Heater

Visual alarm provided

Dual digital temperature display allows easy visualization of patient temperature and set point.

Baby Bed : The spacious baby bed incorporates easy drop down/ fold down glasses for maximum access to the patient.

Castors : The system is fitted with 100mm high quality castors having brakes in front 2 castors.

Additional Features

If needed the unit can be programmed by the user itself for any of the following features:

(a) Count Up / down Timer.

(b) Check baby Alarm.

(c) Pre warming mode.

Dimensions(Approx.)

Technical Specification Electrical Specification

Height- 1850mm Power supply - 200V / 50Hz

Length- 1050mm Heater Power - 650 Watts

Width- 650mm Fuse - 5 Amp.

We are Manufacturer, Exporters and Suppliers of Infant Radiant Warmer Diamond with Phototherapy at wholesale competitive prices from India.

MODEL NO:

BCRW 121 INFANT RADIANT WARMER DIAMOND WITH UNDER SURFACE PHOTOTHERAPY CFL

BCRW 122 INFANT RADIANT WARMER DIAMOND WITH UNDER SURFACE PHOTOTHERAPY LED