

Infant Radiant Warmer Gold without Phototherapy manufactured by Desco Medical India.



Specification:-

Infant Radiant Warmer Gold Series has Advance micro processor controlled system with self check on the electronic circuitry skin / Air Manual control.

The Manual / Servo control mode provides more flexibility while optimizing the thermal environment.

It has Digital display of baby temperature .

In addition to the skin / Air temperature, display is larger than the set temperature display.

The spacious baby bed incorporates easy drop down sides for maximum access to the baby.

The control panel has integrated alarm system with both audible and visual indicators in the event of

• Probe fail • Skin High • Skin Low • Power Fail • Heater Fail

Skin temperature can be converted from Celsius to Fahrenheit .

Baby Bed : The spacious baby bed incorporates easy drop down/ fold down glasses for maximum access to the patient.

Manual Mode : We provide Facility of converting the warmer from servo mode to manual mode for non-stop working of the unit.



Infant Radiant Warmer Gold Without Phototherapy (MODEL NO : BCRW 116)



Castors: The system is fitted with 100mm high quality castors having brakes in front 2 castors. Antistatic castors.

The unit is designed in such a way that the main control system can be swivelled to the other side for access from all the four sides to the patient in critical condition.

Heater Assembly : The heater assembly located above the centre of the support structure consists of an integral examination light and a 650 watts heating system (Silica or quartz / Ceramic Heater) for quick warming and even distribution of heat on the baby bed.

The heater can be swivelled in either direction for taking X-ray and doing procedures.

Option : Weighing scale

Electrical Specification (IEC complaint)

Power supply - 200V / 50Hz

Heater Power - 650 Watts

Fuse - 5 Amp.

Technical Specification

Height 1850mm , Length 1050mm , Width 650mm

We are Manufacturer, Exporters and Suppliers of Infant Radiant Warmer Gold without Phototherapy at wholesale competitive prices from India

MODEL NO:

BCRW 116 INFANT RADIANT WARMER GOLD WITHOUT PHOTOTHERAPY