

Infant Trolley With Resuscitator Manufactured by Desco Medical India



Specification:-

Infant Trolley With Resuscitator is mainly used in delivery room, baby wards, NICU.

Infant Trolley comes with a big shelf at the bottom for keeping consumables.

It can adjust the oxygen concentration from 21% to 100% .

The MAX-P and PIP are safe, controlled and accurate.

PEEP is constant and accurate.

Infant Trolley With Resuscitator is simple to operate and flexible.

The output Pressure would not be influenced by the operator experience, training, attention and fatigue.

Technical Parameter :

Pressure Gauge - 2~10 Kpa (-20~100cm H₂O)

MAX-P- 0~15 LPM 6 Kpa (60cm H₂O)

PIP- 5~15 LPM 0.2~5.7 Kpa (2~57cm H₂O)



Infant Trolley With Resuscitator (MODEL NO : BCRT 101)



The non motorized slow suction can be used with a simple oxygen flow through the oxygen humidifier when connected to an oxygen cylinder / wall unit.

Infant Trolley is manufactured with heavy round pipes.

Acrylic transparent Seamless Crib with mattress.

Body section is mounted on 7.5 cm castors.

It has IV Stand rod.

Oxygen Cylinder Trolley with Humidifier Bottle.

Product Size: 93 x 49 x 80 cm.

Knock Down Construction to save freight.

We are Manufacturer, Exporters and Suppliers of Infant Trolley With Resuscitator at wholesale competitive prices from India.

MODEL NO:

BCRT 101 INFANT TROLLEY WITH RESUSCITATOR