

Estd. 1966



DELIVERING TO THE WORLD

Infusion Blood Warmer (MODEL NO : BBFW 101)



Infusion Blood Warmer Manufactured by Desco Medical India

Specification:-

- * Double hardware and software triple temperature control protection device, ensure warming safety.
- * High temperature forced power-off protection.
- * Warmer speed faster, that is just 2-mins fast heating to the set temperature, save more time, get the precious time for effective blood transfusion;
- * The deviation of heating board within $\pm 1^{\circ}\text{C}$, which is more accurate, safe and heating effectly.
- * The door buckle design can prevent the door from accidentally opening during the heating process and affect the accuracy of heating, add the alarm when the door open.
- * Various temperature control protection technology, $^{\circ}\text{C}$ and F unit convert, convenient for usage habits in different regions.
- * Two fixing methods of card holder and sling.
- * There is silent alarm function, which ensure quiet treatment environment.
- * Design small, making exquisitely, as light 0.2kg.
- * Brightness display, real-time temperature can be watched from a distance can.



Infusion Blood Warmer (MODEL NO : BBFW 101)



Technical Specification :

Temperature Setting	°C: 35-42°C °F: 95.0°F-107.6°F
Temperature Accuracy	°C: ±1°C °F: ±1.8°F
Increment	°C: 0.1°C °F: 0.2°F
Optional Unit	°C/°F
Warming up time	≤2min
Overheat Protection	46°C
Alarm	Over heat, low temperature, system error, door open, warming up alarm
Display	LED
Power Supply	Input: 100-240V(AC power), 50Hz/60Hz, 0.3-0.7A Output: 12V(DC Power) 2.5A"
Power Consumption	85VA
Operation Environment	Temperature: +5°C-+30°C Relative humidity: 20%-80% Atmospheric pressure: 70.0Kpa-106.0Kpa
Transportation and Storage Environment	Temperature: -20°C-+55°C Relative humidity: 10%-93% Atmospheric pressure: 50.0Kpa-106.0Kpa
Degree of protection against ingress of liquids	Class II, IPX1
Dimension	Size: 180mm*70mm*35mm
Weight	Without Knob: 0.2KG With Knob: 0.4KG"

We are Manufacturer, Exporters and Suppliers of Infusion Blood Warmer at wholesale competitive prices from India.

Model No

BBFW 101

Infusion Blood Warmer