

Estd. 1966



DELIVERING TO THE WORLD

Ophthalmic OT Table Manual (MODEL NO : OTOP 101)



Ophthalmic OT Table Manual Manufactured by Desco Medical India



Specification:-

Three sectional table is provided with drainage system in the head section. The table has sturdy base, column made up of stainless steel, antibacterial and is easy to clean. Electromotive height and trendelenburg positions are with the help of light weight handheld control unit. Easy detachable leg and head section.

Features :

- * The base is eccentric.
- * Stainless steel fitting.
- * Three sectional Radio-translucent top.
- * Base and cover are made of attractive , easy to clean stainless steel.
- * Sophisticated mechanics provide smooth, stepless articulation of the table top for precise patient positioning.
- * Leg & Head section are interchangeable.



Ophthalmic OT Table Manual (MODEL NO : OTOP 101)



Technical Specs:-

Length of Table	1900 mm
Width of Table Top	520 mm
Minimum Height (without Mattress)	Min. 750 mm
Maximum Height (without Mattress)	Max. 1000 mm
Trendelenburg	25°
Reverse Trendelenburg	25°
Leg Section	0°-90°

We are Manufacturer, Exporters and Suppliers of Ophthalmic OT Table Manual at wholesale competitive prices from India.

MODEL NO:

OTOP 101 Ophthalmic OT Table Manual