

## Suction Catheter Closed System Manufactured by Desco Medical India



### Specification:-

- Suction Catheter Closed System is made from non-toxic, medical grade PVC.
  - Suction Catheter Closed System is used to maintain ventilation and oxygen therapy while suctioning in mechanically ventilated patients.
  - Minimises the risk of Ventilator Associated Pneumonia (VAP).
  - Kit comes with 3 easily replaceable catheters, minimizing the need to replace entire suction system every 24-72 hrs.
  - Catheter with lateral eyes, smooth tip and depth markings for maximum patient comfort.
  - MDI port provides an effective, quick and convenient way of drug delivery.
  - Irrigation port with one way valve allows lavage and irrigation of the suction catheter.
  - On/off valve allows suction catheter tip to be isolated and creates an enclosed cleaning chamber.
  - Transparent connectors and durable sleeve facilitates catheter manipulation and prevents infiltration of blood and sputum.
  - 360° swivel elbow enables rotation for enhanced patient comfort and reduces torque.
- We are Manufacturer, Exporters and Suppliers of Suction Catheter Closed System at wholesale competitive prices from India.

Estd. 1966



DELIVERING TO THE WORLD

# Suction Catheter Closed System

## (MODEL NO : MDDC 201)



<b>Model No</b>	<b>Sizes in FR</b>	<b>Outer Dia. In MM</b>	<b>Inner Dia. In MM</b>	<b>Length in MM</b>	<b>Colour Code</b>
MDDC 201	6	2.0	1.4	300	Pea Green
MDDC 202	8	2.7	1.8	300	Cambridge Blue
MDDC 203	10	3.3	2.4	600	Black
MDDC 204	12	4.0	2.8	600	White
MDDC 205	14	4.7	3.2	600	Green
MDDC 206	16	5.3	3.8	600	Orange