

## Traction Unit Computerised with Bed Supplied by Desco Medical India

  
**CE**

  
**D&B**  
Decide with Confidence

Estd. 1966

  
**DESCO**  
DELIVERING TO THE WORLD

### Traction Unit Computerised with Bed

  
**Cervical Traction**



[www.descomedicalindia.com](http://www.descomedicalindia.com)

[www.descoonline.in](http://www.descoonline.in)

### Specification:-

Hi-Traction micro computer controlled unit for Intermittent / Constant, Lumbar - cum - Cervical Traction is based on solid- state circuitry.

The three section fixed height powder coated Bed with storage cabin underneath is an integral part of the time.

Technical :

Therapy Modes : EIGHT Traction Modes for either Lumbar or Cervial Traction, in progressive, Regressive, Intermittent & Static Modes etc.

Memory / Program : 30 Treatment Patterns can be programmed & stored in non-volatile memory. 24 Programers are feeded and six programs are Free).

Traction Force : Adjustable from 05kg to 60kg.

Base Force : Adjustable from 00kg to 54kg.



# Traction Unit Computerised with Bed (MODEL NO : PTTA 121)



## Traction Unit Computerised with Bed Supplied by Desco Medical India

### Specification:-

Auto-Tension Adjustment : The Set Traction force is automatically maintained throughout the operations.

Traction Hold Time : Adjustable From 00 to 99 sec.

Base Traction Hold Time : Adjustable From 00 to 99 sec.

Treatment Time : Continuously adjustable from 00 to 99 minutes, with acoustic signal and automatic reduction of traction force.

Displays : Digital Display of Traction Force, Base Force, TF Hold Time, BF Hold Time, Treatment Time & Graphical Displays of Traction Mode.

### Traction / Patient Safety :

1. Automatic System Check in Software on Start of Traction.
2. High Force Signal is Activated if Traction Force is set above 18kg & Acceptance of High-Force is required before start of treatment.
3. Automatic Tension release when emergencies occur. Patient can stop with hand held switch.
4. Display of Error Contents.
5. Power Supply : 220-240V AC. (use Stabilizer)

### Traction Table :

1. The Fixed Height Upholstered traction table is ideal for Horizontal Cervical and Pelvic / Lumbar Traction.
2. The Lumbar Section Rolls freely on nylon rollers for friction free lumbar pull. This section can also be locked in any position if desired. Head section has an Adjustable incline.
3. Cervical Traction angle is controlled by adjustment of the angle adjustment bracket.
4. Size : Table is 212cm X 68cm X 63cm high and comes complete with storage shelf.

### Accessories :

Adjustable Flexion Stool, Head Halter, Pelvic & Thoracic Belts.

We are Exporters and Suppliers of Traction Unit Computerised with Bed at wholesale competitive prices from India.

### MODEL NO:

**PTTA 121**      Traction Unit Computerised with Bed