

## Traction Unit with Fixed Height Table Supplied by Desco Medical India



### Specification:-

Microprocessor based "Nice-Trac" unit with Continuous & Intermittent traction modes is an accurate and sturdy unit with three section fixed height Traction Table.

Traction machine works on solid state circuitry with reliable, silent & smooth geared motor.

### Technical :

Traction Force : Direct Pulling 05 to 45kg in Twelve steps.

Maximum Traction Force : With Doubler Pulley upto 90kgs.

Hold Time : 00 to 99 Seconds, adjustable in 1 second Step.

Rest Time : 00 to 99 Seconds, adjustable in 1 second Step.

Treatment Time : 00 to 60 minutes, adjustable in 1 minute Step.

Alarms : At the end of treatment time.

Safety : Patient safety switch provided with Alarm (Audio & Video).

Power Supply : 220V AC



# Traction Unit with Fixed Height Table (MODEL NO : PTTA 101)



## Traction Unit with Fixed Height Table Supplied by Desco Medical India

### Specification:-

Mounting : Can be mounted on Table or on Wall.

The fixed height upholstered traction table is ideal for Horizontal Cervical and Pelvic / Lumbar Traction.

The lumbar section rolls freely on nylon rollers for friction free lumbar pull. This section can also be locked in any position if desired.

Cervical traction angle is controlled by vertical adjustment of the Machine mounting board.

Size : Table is 182cm long x 66cm wide x 76cm high and comes complete with Storage Shelf.

### ACCESSORIES:

Cervical Head Holder with Bar.

Lumber Traction Belts with Bar.

Adjustable Flexion Stool, Padded.

Main Cord.

Weight Doubler Pulley.

We are Exporters and Suppliers of Traction Unit with Fixed Height Table at wholesale competitive prices from India.

### MODEL NO:

**PTTA 101**      Traction Unit with Fixed Height Table