

Traction Unit with Hi Low Table Supplied by Desco Medical India



Specification:-

An Intermittent traction modes is an accurate and sturdy unit works on solid state circuitry with reliable ,silent & smooth geared motor. Unit comes with with three section Hi-Low Traction Table.

Technical :

- Traction Force : Direct Pulling 05 to 45kg in Twelve steps.
- Maximum Traction Force : With Doubler Pulley upto 90kgs.
- Hold Time : 00 to 99 Seconds , adjustable in 1 second Step.
- Rest Time : 00 to 99 Seconds, adjustable in 1 second Step.
- Treatment Time : 00 to 99 minutes, adjustable in 1 minute Step.
- Alarms : At the end of treatment time.
- Safety : Patient safety switch provided with Alarm (Audio & Video).
- Operating Voltage : 220-240V AC / 50Hz



Traction Unit with Hi Low Table (MODEL NO : PTTA 111)



Traction Unit with Hi Low Table Supplied by Desco Medical India

Specification:-

Hi-Low Traction Table :

Traction Table : Hi-Low Electric operated, Three section table is suitable for traction / treatment and manipulation.

Height Adjustment : Electrically with Hand-Switch from 50cm to 95cm

Upholstered Top Size : 195cm Long X 70cm Wide

Traction Mount Platform : Adjustable in Height for proper pull.

Sections : Three Sections, drop end section with breathing hole and plug adjusts from -20 to +20 deg. Non-Friction rolling Lumbar section Tilts 60 degree above horizontal.

Traction Support : Traction Unit can be bolted on mounted board, for cervical and lumbar traction.

ACCESSORIES:

Cervical Head Holder with Bar.

Lumber Traction Belts with Bar.

Adjustable Flexion Stool, Padded.

Main Cord.

Weight Doubler Pulley.

We are Exporters and Suppliers of Traction Unit with Hi Low Table at wholesale competitive prices from India.

MODEL NO:

PTTA 111 Traction Unit with Hi Low Table