

## Un-Weigh Mobility Trainer with Treadmill Supplied by Desco Medical India

  
**CE**

  
**D&B**  
Decide with Confidence

Estd. 1966

  
**DESCO**<sup>®</sup>  
DELIVERING TO THE WORLD

### Un-Weigh Mobility Trainer with Treadmill



The image shows a white and blue Un-Weigh Mobility Trainer with a treadmill. It features a central vertical column with a digital control panel and a suspension bar at the top. A black harness is attached to the suspension bar. The treadmill is mounted on a base with four casters for mobility.

[www.descomedicalindia.com](http://www.descomedicalindia.com)

[www.descoonline.in](http://www.descoonline.in)

### Specification:-

The Un-weighing Trainer enables partial weight - bearing therapy to be conducted with the assurance of patient Comfort and Safety, and with convenient access to the patient for manual observation and assistance.

The electrical un-weighing trainer is designed to apply vertical support to remove the stress of bearing body weight, unit digitally controls weight bearing, and promotes proper posture and balance over a treadmill or the ground.

Harnesses give security to patients with limited trunk strength. The adjustable suspension bar adjusts to accommodate tall, Children & Adults.

## Un-Weigh Mobility Trainer with Treadmill Supplied by Desco Medical India

### Specification:-

#### Technical :

- \* Un-weighing : Can be set 0% to 100% of the body weight, digitally.
- \* Un-weighing System : Electrical the unit comes with mains operated built in lift with battery backup.
- \* Vertical Oscillation Upto 24", ideal for sit to stand exercises also.
- \* Suspension : Single point with 360 degree rotational bar to help patient turn easily.
- \* Control Panel : Digitally display, Patient body weight and Un-weighing weight.
- \* Capacity : 0-120 kg un-weighing range.
- \* Un-Weighing Harness : Two Small and Standard size.
- \* Handrails : For patient's support, Adjustable & Removable.
- \* Frame measures : 42"L x33" wide inside & fitted with Four ,6" lockable wheels.
- \* Power Supply : Unit Operates on 220VAC, with battery backup system.

We are Exporters and Suppliers of Un-Weigh Mobility Trainer with Treadmill at wholesale competitive prices from India.

#### MODEL NO:

**PTTF 161**

Un-Weigh Mobility Trainer with Treadmill