



# Un-Weigh Mobility Trainer without Treadmill (MODEL NO : PTTF 151)



## Un-Weigh Mobility Trainer without Treadmill Supplied by Desco Medical India



Decide with Confidence

### Un-Weigh Mobility Trainer without Treadmill



[www.descomedicalindia.com](http://www.descomedicalindia.com)

[www.descoonline.in](http://www.descoonline.in)

#### Specification:-

The Un-weighing Trainer enables partial weight - bearing therapy to be conducted with the assurance of patient Comfort and Safety, and with convenient access to the patient for manual observation and assistance.

The electrical un-weighing trainer is designed to apply vertical support to remove the stress of bearing body weight, unit digitally controls weight bearing, and promotes proper posture and balance over a treadmill or the ground.

The system allows patient to stand upright and use both arms freely. Harnesses give security to patients with limited trunk strength.

The adjustable suspension bar adjusts to accommodate tall, Children & Adults.



# Un-Weigh Mobility Trainer without Treadmill (MODEL NO : PTTF 151)



## Un-Weigh Mobility Trainer without Treadmill Supplied by Desco Medical India

### Specification:-

#### Technical :

- \* Un-weighing : Can be set 0% to 100% of the body weight, digitally.
  - \* Un-weighing System : Electrical the unit comes with mains operated built in lift system.
  - \* Vertical Oscillation Upto 24", ideal for sit to stand exercises.
  - \* Suspension : Single point with 360 degree rotational bar to help patient turn easily.
  - \* Patient body weight and Un-weighing weight.
  - \* Control Panel : Digitally display
  - \* Capacity : 0-120 kg un-weighing range.
  - \* Un-Weighing Harness : Two Small and Standard size.
  - \* Handrails : For patient's support, Adjustable & Removable.
  - \* Therapist Seat : For therapist to sit during assistive exercises.
  - \* Frame measures : 42"L x 33" wide inside to accommodate Treadmill & fitted with Four ,6" lockable wheels.
  - \* Power Supply : Unit Operates on 220VAC, with battery backup system.
  - \* Note : Assembled working size of unit ; Height 92" x Width 39" x Length 43" approx.
- We are Exporters and Suppliers of Un-Weigh Mobility Trainer without Treadmill at wholesale competitive prices from India.

### MODEL NO:

**PTTF 151**

Un-Weigh Mobility Trainer without Treadmill