

Unolock Syringe with Needle (MODEL NO : MDSN 501)



Dispovan Syringe with Needle Manufactured by Desco Medical India



Specification:-

- Luer lock Single use syringes with needle.
- Secure fitment of needle due to presence of luer lock.
- High clarity plunger and barrel.
- Better plunger retention due to locking at the base.
- Barrels are made of nontoxic, medical-grade polypropylene compatible with any medication.
- The tip of the barrel is luer lock having a 6% luer taper.
- Plungers are made of nontoxic, medical-grade polypropylene compatible with any medication.
- Gaskets are made of natural rubber which is chemically inert and compatible with medication for short term contact drug delivery application.
- Luer Lock nozzle present in Unolok syringes ensures secure fitment of the needle.
- We are Exporters and Suppliers of Unolock Syringe with Needle at wholesale competitive prices from India.



Unolock Syringe with Needle (MODEL NO : MDSN 501)



Model No	Capacity	Gauge	Length in mm	Dia in mm	Type of Nozzle
MDSN 501	2 ml	21G	25 mm	0.80	Luer Mount with Central Nozzle
MDSN 502	2 ml	21G	38 mm	0.80	Luer Mount with Central Nozzle
MDSN 503	2 ml	22G	25 mm	0.70	Luer Mount with Central Nozzle
MDSN 504	2 ml	22G	32 mm	0.70	Luer Mount with Central Nozzle
MDSN 505	2 ml	22G	38 mm	0.70	Luer Mount with Central Nozzle
MDSN 506	2 ml	23G	25 mm	0.60	Luer Mount with Central Nozzle
MDSN 507	2 ml	23G	38 mm	0.60	Luer Mount with Central Nozzle
MDSN 508	2 ml	24G	25 mm	0.55	Luer Mount with Central Nozzle
MDSN 509	3 ml	21G	25 mm	0.80	Luer Mount with Central Nozzle
MDSN 510	3 ml	21G	38 mm	0.80	Luer Mount with Central Nozzle
MDSN 511	3 ml	22G	25 mm	0.70	Luer Mount with Central Nozzle
MDSN 512	3 ml	22G	32 mm	0.70	Luer Mount with Central Nozzle
MDSN 513	3 ml	22G	38 mm	0.70	Luer Mount with Central Nozzle
MDSN 514	3 ml	23G	25 mm	0.60	Luer Mount with Central Nozzle
MDSN 515	3 ml	23G	38 mm	0.60	Luer Mount with Central Nozzle
MDSN 516	3 ml	24G	25 mm	0.55	Luer Mount with Central Nozzle
MDSN 517	5 ml	21G	25 mm	0.80	Luer Mount with Central Nozzle
MDSN 518	5 ml	21G	38 mm	0.80	Luer Mount with Central Nozzle
MDSN 519	5 ml	22G	25 mm	0.70	Luer Mount with Central Nozzle
MDSN 520	5 ml	22G	32 mm	0.70	Luer Mount with Central Nozzle
MDSN 521	5 ml	22G	38 mm	0.70	Luer Mount with Central Nozzle
MDSN 522	5 ml	23G	25 mm	0.60	Luer Mount with Central Nozzle
MDSN 523	5 ml	23G	38 mm	0.60	Luer Mount with Central Nozzle
MDSN 524	5 ml	24G	25 mm	0.55	Luer Mount with Central Nozzle
MDSN 525	10 ml	20G	38 mm	0.90	Luer Mount with Central Nozzle
MDSN 526	10 ml	21G	38 mm	0.80	Luer Mount with Central Nozzle
MDSN 527	10 ml	22G	32 mm	0.70	Luer Mount with Central Nozzle
MDSN 528	10 ml	22G	38 mm	0.70	Luer Mount with Central Nozzle
MDSN 529	20 ml	20G	38 mm	0.90	Luer Mount with Central Nozzle
MDSN 530	20 ml	21G	38 mm	0.80	Luer Mount with Central Nozzle
MDSN 531	20 ml	22G	32 mm	0.70	Luer Mount with Central Nozzle
MDSN 532	20 ml	22G	38 mm	0.70	Luer Mount with Central Nozzle