

Estd. 1966



DELIVERING TO THE WORLD

# Volumetric Infusion Pump (MODEL NO : IEIP 101 )



## Volumetric Infusion Pump Supplied by Desco Medical India



Decide with Confidence

### Volumetric Infusion Pump

Estd. 1966



DELIVERING TO THE WORLD



[www.descomedicalindia.com](http://www.descomedicalindia.com)

[www.descoonline.in](http://www.descoonline.in)

Estd. 1966



DELIVERING TO THE WORLD

# Volumetric Infusion Pump (MODEL NO : IEIP 101 )



## Specification:-

Infusion pumps may be capable of delivering fluids in large or small amounts, and may be used to deliver nutrients or medications – such as insulin or other hormones, antibiotics, chemotherapy drugs, and pain relievers.

The infusion pump has audible and visual alarm for occlusion, empty, low battery, end of infusion, door open, wrong setting etc ,

HD LCD Display, high capacity words, friendly user interface, dynamically display work status .

The infusion pump is compatible with any brand of infusion sets after correct calibration ,

The infusion pump can preset solution volume, which can greatly reduce the workload of nurses

The infusion pump has two work modes: ml/h and drop/min can switch freely

There are three levels of occlusion: high, middle and low for the pump

The pump has purge function

The infusion pump will automatically open up KVO as infusion is completed.

This infusion pump also can automatically record the settings for last transfusion

Infusion Rate : 1ml/h - 1200ml/h

Purge Rate : 100ml/h - 1200ml/h

Volume Limit : 0ml- 9999ml

Speed Range ml/h-1200ml/h (increment of 1ml/h) .

Infusion Accuracy with in  $\pm 5\%$

High Speed Infusion 100ml/h--1200ml/h (can be calibrated) ,

Volume Limit 1ml-999.9ml roh increment of 0.1ml

Infusion Volume Done 9999ml(increment of 0.1ml)

Weight:-1.8kg

Dimension:-140 mm (L) x 157mm (W) x 220 mm (H)

We are Exporters and Suppliers of Volumetric Infusion Pump at wholesale competitive prices from India.

**MODEL NO**

**IEIP 101**

**Volumetric Infusion Pump**